


Combination Therapy OASYS with CPAP Pro

In some cases, Oral Appliance Therapy alone will not reduce the AHI satisfactorily for the Sleep Physician. Combination Therapy can be incorporated, by modifying the appliance with a removable bracket, which will allow the patient to wear the OASYS appliance with CPAP or use the appliance as originally designed, for travelling or when ill. Using the OASYS appliance with the CPAP Pro, minimally invasive Nasal Pillow Assembly, can also reduce CPAP pressure requirements and provide more comfort than other masks.

Fixed for Home Use Removable for Travel



Removable Bracket



A custom bracket has been designed to attach the OASYS to the CPAP Pro provided. Hooks and Elastics are included to encourage mouth closure. A Sleep Physician or Sleep Technician will need to determine the appropriate air pressure required, and possibly modify the CPAP Pro for comfort and proper function.



720 Sunrise Avenue, #230-A * Roseville, CA 95661
(916) 865-4528 www.dreamsystemsdentallab.com